## **Photon Sound Beam compared to Other Instruments:**

## **PSB vs Rife Technology**

### PSB vs Photon Genie, Photon Genius, LBG, and the BCX Ultra

## **Photon Sound Beam vs Rife Technology**

Since 1988 we have been educating, developing, and promoting a wide variety of Vibrational Energy Instruments. We have extensive background in both allopathic and holistic medicine and offer full health consultations and technical support. Our focus has always been to evolve an instrument beyond basic Rife Technology. This is why we have been continually advancing the instrumentation of the *Photon Sound Beam* -- now a time-tested comprehensive Lymphatic Rife & Lakhovsky Technology.

In Electro-Medicine *delivery* is everything! The key is to use *direct-contact applicators* with biological plasma-electricity AND high-harmonic frequencies to fully penetrate the body's subtle and gross energy pathways which enhances the delivery and bio-availability of a treatment.

The internet is full of information on Rife Technology created by many different makers, but they do not use *Electric Plasma Tubes* to deliver Rife Frequencies and other frequencies -- which is the core value of the *Photon Sound Beam (PSB)*. "Rife Technology" as most people understand does not go beyond the allopathic model to simply eliminate the pathogens. The *Photon Sound Beam* is allopathic *and* holistic. It eliminates the pathogens *while* deeply cleansing and nourishing the body for long-term regenerative & rejuvenating values.

#### There are 2 types of Rife Technology:

- Classical Rife/Bare Broadcast Instruments -- modern versions of Rife's Beam Tube which use Radiowave Plasma Tubes to broadcast at a distance
- Rife/Crane Frequency Generators which use skin-contact Metal Rods or Metal/Gel Pads to conduct an electrical current which does not match the extremely low micro-current of the body

Classical Rife/Bare Broadcast: There are severe limitations with any broadcast instrument with signal loss and signal distortion through our electromagnetic & cell phone-polluted environment -- even with a CB rectilinear amplifier. Also, this high wattage radiowave broadcast may be harmful to delicate heat-sensitive neurology. We and several others including a kinesiological group of 30 years in San Diego have always tested a deleterious effect with a body response to Rife/Bare broadcast instruments.

Rife/Crane Frequency Generators: We have referred to this Rife Technology as "dinosaur technology" since the early 1990's because of a frequency impedance problem.\* This Metal Rod/Pad delivery system uses an electrical current which does not match that of our bodies and it can only travel along the main energy channels without adequate penetration and without any lymphatic cleansing.

\*Ed Skilling and Robert Beck first made contemporary Metal Rod/Pad units in the 1980's and the specs were published for anyone to copy. When we met with Ed Skilling in the late 80's he said that metal (or gel pad) to-skin-contact "distorts the frequency because of a skin impedance problem". That's why he introduced his close-contact *RF Sound Probe* in 1988 which is now a component of the *Photon Sound Beam*.

The primary difference between basic Rife Technology and the *Photon Sound Beam* is that these 2 types of "Rife Technology" do not provides *lymphatic & tissue cleansing* which is fundamental to *any* health issue. Secondly, it has 2 components which work in synergy to deliver extremely pure and bio-available Rife Frequencies, brainwave and other frequencies. It expands an aggressive allopathic technology to include cleansing and nourishing for more universal, curative, preventative and creative health benefits.

#### The new *Photon Sound Beam* is a comprehensive Lymphatic Rife and Lakhovsky Instrument.

The *Photon Sound Beam* effectively delivers a full-harmonic spectrum of Rife, brainwave and other frequencies with 2 Components:

- Noble Gas Plasma Tubes (Electrically Ionized & RF Enhanced, Xenon/Argon/Krypton)
- Radio Frequency Sound Probe (Dual Crystal Oscillators, Lakhovsky Technology)

NOTE: The **PSB XII** is our most popular personal-use /clinical model. The **PSB Infinity RF** is our clinical model, and adds a Radiowave Ionization option with full-frequency selection and programming.

The *Photon Sound Beam* delivers *plasma-electricity* and *radiowave signals* directly to every cell in the body with direct skin-contact glass applicators. This quickly removes all the blockage & barrier issues and opens-up all the energy and information pathways. It deeply cleanses and oxygenates the body's tissues & fluids while de-congealing stagnant lymph. And this allows the frequencies to effectively integrate into the body to rapidly affect change. This also prevents a negative detoxification response from pathogen die-off during a treatment. Plasma-Electricity delivers electricity at the level of the cell by a process called electrical disassociation which restores the electrically negative (-) charge to cell membranes which is weakened from acid conditions in the body.

At the same time the *PSB's RF Sound Probe* also devitalizes the pathogens and re-establishes the ideal cellular resonance and polarity of all the *healthy cells*. This unique component of the *Photon Sound Beam* is based upon the lengthy scientific research of Georges Lakhovsky in the 1920's, who demonstrated -- well before Raymond Rife's work -- that a radiowave frequency broadcast with *high harmonics* always proved better in treating cancer. The *RF Sound Probe* delivers a harmonically-rich radiowave broadcast applied *directly* to the body or at a close distance. The *RF Sound Probe* has 2 quartz crystal oscillators (in Phi Ratio) and delivers extremely pure & high-harmonic frequencies deep into the body with a unique, highly-scalar RF broadcast antenna.

Our completely upgraded *Photon Sound Beam* is a comprehensive wrap-around technology and are a profound improvement over any type of Rife Technology. This is why since 1994 our vastly upgraded *Photon Sound Beam* has consistently provided a rapid and universal response from anyone, for any issue, at any level of health!

# Photon Sound Beam vs Other Instruments

#### Photon Genie, Genius, TeslaStar, LBG, BCX Ultra

The *Photon Genie* was created in the late 1990's and is essentially a re-make of the prototype 1994 *Inert Gas Sound Probe* designed by Ed Skilling in 1994. At this same time it was re-named the *Photon Sound Beam* by *Integrated Healing Systems* having promoting Ed Skilling's *Sound Probe* beginning in 1988. The *Sound Probe* is now an upgraded *component* of all our *Photon Sound Beam* models.

In the late 1990's, Ed Skilling hired an engineer to add an additional output coil for the glass Plasma Tubes to a 1994 PSB Prototype unit and called it the Photon Genie. By this time we had upgraded our technical base with a team of engineers from Sedona Arizona -- who have been dramatically evolving the 1994 Prototype PSB ever since. We now offer a full line of plasma-electrical/RF Photon Sound Beams. (See Photon Sound Beams: Contemporary History)

The *Photon Genie* has an additional high voltage coil but otherwise it has the same circuitry as the *1994 PSB Prototype* still be sold by other companies and still absent basic control features common to most contemporary electro-medical instruments.

The *Photon Genie* is factory-manufactured in Japan and its users have reported no support for repairs, and consequently have had to buy a new unit at price exceeding \$3000. *Integrated Healing Systems* offers lifetime technical support and a 2 year warranty. Our most popular (clinical/personal-use) *Photon Sound Beam XII* has much more evolved instrumentation for *100's of dollars less*.

Since the early 1990's the Photon Sound Beam has been continually evolved now with over 20 upgrades in circuitry, control features and hardware, including it heavy duty, triple gas Bi-Phasic Plasma Tubes. It has been greatly advanced in overall value and long term benefits by incorporating important control features such as Pulse Rate Control, Frequency Selection and Frequency Sweep which are important in enhancing regeneration and most importantly in overcoming body accommodation -- a fundamental problem in electro-medicine.

Our personal-use/clinical **Photon Sound Beam XII** our clinical **PSB Infinity RF** have been our premier instruments since 2005 and were preceded by our **PSB VIII**, **PSB VII** and the **1994 PSB Prototype**.

#### The Photon Genie, Photon Genius, and the Light Beam Generator have:

- No Frequency Selection (728 Hz Universal Rife Frequency ONLY)
- No Pulse-rate Control (found to enhance regeneration)
- Argon-ONLY Plasma Tubes the PSB has a Triple Gas Tube (Xenon/Argon/Krypton)
- No advancements in Noble Gas Applicators -- Both the LBG and Photon Genie are limited by the use off-the-shelf glass tubes. The LBG output has limited penetration (wattage) of the cup- shaped glass applicator. NOTE: The Genie uses painted glass tubes for use in neon signs which obstructs the transmission of the elemental color frequencies and subtle energies generated from ionizing the Noble Gases. The "magic" of this type of any Plasma Electrical instrument is the use of the remarkable media inside the glass applicators for transmitting frequencies: a near vacuum-state with electricity and Noble Gases.

ABSENT in the *Photon Genie, Photon* Genius, *BCX Ultra* and the *Light Beam Generator* — The *Photon Sound Beam XII* has:

- Electrically Ionized Glass Plasma Tubes with 3 Noble Gases (Xenon/Argon/Krypton) for lymphatic and tissue cleansing
- 12 Frequency Selections including a 36 Frequency Sweep of Priority Rife & Brainwave Frequencies, and dial selections for 4 Universal Rife and 5 Classical Cancer/Viral Frequencies, Alpha and Theta
- Advanced *RF Sound Probe* (called the *Flat Pack* for the *Genie*) -- the most harmonically-rich frequency emitter of any technology available, in our opinion. It has 2 quartz crystal oscillators at 25 and 40 Megahertz -- broadcasting 2 *Carrier RF signals* in Phi Ratio which carry harmonically-rich Sound Frequencies, Brain Wave Frequencies, and Rife Frequencies up to 10,000 Hz. *It also uses an advanced Scalar Broadcasting Antenna*.
- Individual Pause and Pulse Controls can be used to harmonize with breath and heart rates and prevents body-accommodation, Pulsation Rates have been shown to affect regeneration rates in experiments with salamander neurology and the re-growing of tails.
- Heavy Duty construction including non-conductive components such as hardwood console sides to maximize plasma-electrical output, and extra thick high voltage cords & protectors
- **20 Upgrades in Circuitry and Features** since 1996 beyond 94 *PSB Prototype* & *Photon Genie*

The **Photon Genius** is a wall broadcast plasma-tube instrument, and it has the same limitations as the **Photon Genie** by the same manufacturer. Like the **Genie**, there is only a single frequency with no expanded control features. Many users report the direct-contact **Photon Genie** to be a more effective instrument.

NOTE: ANY broadcast instrument such as the Photon Genius has severe problems with signal-loss and signal distortion even across a few feet distance, especially in our cell phone/EMF- polluted world. With broadcasting, independent testing has shown that lymphatic drainage is severely minimized compared to direct contact Plasma Tubes. When you double the broadcast distance the signal-loss issquared.

These same broadcast limitations distortion also apply to the *Vibe Machine, Novalite* and even the upgraded version *TeslaStar*. These instruments are stand-up broadcasters which have large diameter plasma-tubes and a Pulsed EMF coil. They are nice additions in a clinical setting, but they are up to 4 times the cost AND again, *direct contact glass plasma-tubes* always deliver frequencies and cleanse the fluids & tissues much more purely and effectively. The Pulsed Magnetic Field component in general provides more of a "spa effect" rather than quickly addressing symptoms. Magnetic field broadcasting is *severely* weakened with distance.

The **Light Beam Generator** and the high-priced *ST-8* model lack sufficient focused penetration with an array of lower wattage plasma-tube applicators. It has a more diffuse effect on the lymphatic system AND lacks the *RF Sound Probe* component,

The **BCX Ultra** uses radio frequencies to ionize the Plasma Tubes which does not provide the fundamental value of the *Photon Sound Beam* and ANY effective vibrational energy technology which is *Plasma-Electricity* from Electric Plasma Tubes NOTE: *RF ionization* is a component

already included in the *PSB Infinity RF*. A second *BCX Ultra* component is a Metal Rod/Pad Rife/Crane instrument which does not clean the lymph and tissues and distorts the frequencies. *See <u>PSB vs Rife Technology</u>*. We have added this component to the *PSB Infinity RF* with our special order *Lifepulse Pro* for those insisting on including it.

Our clinical model, the **PSB Infinity RF** adds additional components to the **PSB XII.** It adds Radiowave Ionization of the glass Plasma-Tubes, an LCD display, and full-frequency selection & programming, and includes the <u>Electroherbalism Frequency Lists</u> book. The unique RF Ionization of the Plasma-Tubes allopathically (more aggressively) eliminates the parasites and pathogens with deep frequency penetration. The Infinity RF is highly recommended for chronic or acute issues or for practitioners treating a variety of conditions. Users have reported the extra benefit of enhanced RF-penetration of frequencies using 60 Selections of Frequency Programs especially in urban EMF & frequency -polluted settings.

Integrated Healing Systems continues to offer lifetime technical support of our technologies, as well as a 2 year warranty, and holistic health consultations. We are always available to discuss any technology.

We are fortunate to have had contact with practitioners and clients from all over the world during these last few decades. This clinical feedback, experiential knowledge and the contribution by our amazing engineers has generated an effective plasma-electrical Rife/Lakhovsky technology -- instruments which are which are easy-to-use and very fast-acting on a wide range of issues for a reasonable price.